


Product Datasheet

beta-Catenin (Phospho-Ser33/37/Thr41) antibody (orb336576)

Description	Rabbit polyclonal antibody to beta-Catenin (Phospho-Ser33/37/Thr41)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3020
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...