


## Product Datasheet

### Aquaporin 7 antibody (orb336776)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to Aquaporin 7 |
| <b>Species/Host</b>        | Rabbit                                    |
| <b>Reactivity</b>          | Human                                     |
| <b>Conjugation</b>         | Unconjugated                              |
| <b>Tested Applications</b> | WB  |
| <b>Note</b>                | For research use only                     |
| <b>Clonality</b>           | Polyclonal                                |
| <b>Dilution Range</b>      | WB: 1:1000-3000                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.           |



Western blot analysis of Jurkat whole ce...