


## Product Datasheet

### GAPDHS antibody (orb338331)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GAPDHS |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human                                |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | WB                                   |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Dilution Range</b>      | WB: 1:1000-3000                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |



Western blot analysis of HeLa whole cell...