

---

## Product Datasheet

### Focal Adhesion Kinase (phospho-Y407) antibody (orb338933)

**Description**

Rabbit polyclonal antibody to PTK2

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Primate, Rat, Zebrafish

**Conjugation**

Unconjugated

**Tested Applications**

IH, WB

**Immunogen**

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y407 of human Focal Adhesion Kinase protein. The exact sequence is proprietary.

**Target**

PTK2

**Preservatives**

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

**Form/Appearance**

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Clonality**

Polyclonal

**Source**


Rabbit

**Uniprot ID**
[P34152](#), [O35346](#), [Q05397](#)
**Entrez**
[5747](#), [14083](#), [25614](#)
**Dilution Range**


WB: 1:500-1000, IHC-P: 1:100-200

**Expiration Date**

12 months from date of receipt.



Western blot analysis of Focal Adhesion ...



Immunohistochemical analysis of Focal Ad...