

## Product Datasheet

### FOXO1 (phospho-S256) peptide (orb339213)

<b>Description</b>	The peptide is used to FOXO1 (pS256) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Porcine, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	FOXO1
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338819 (FOXO1 (phospho-S256) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q12778</a> , <a href="#">G3V7R4</a> , <a href="#">Q9R1E0</a>
<b>Entrez</b>	<a href="#">2308</a> , <a href="#">84482</a> , <a href="#">56458</a>
<b>Expiration Date</b>	12 months from date of receipt.