

Product Datasheet

Estrogen Receptor alpha (phospho-S104) peptide (orb339218)

Description	The peptide is used to Estrogen Receptor alpha (pS104) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ESR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338824 (Estrogen Receptor alpha (phospho-S104) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19785 , P03372 , P06211
Entrez	2099 , 24890 , 13982
Expiration Date	12 months from date of receipt.