

Product Datasheet

VEGFR2 (phospho-Y1214) peptide (orb339221)

Description	The peptide is used to VEGFR2 (pY1214) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	KDR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338827 (VEGFR2 (phospho-Y1214) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35918 , O08775 , P35968
Entrez	25589 , 3791
Expiration Date	12 months from date of receipt.