

Product Datasheet

SHPTP1 (phospho-Y536) peptide (orb339222)

Description	The peptide is used to SHPTP1 (pY536) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PTPN6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338828 (SHPTP1 (phospho-Y536) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29350 , P81718 , P29351
Entrez	116689 , 15170 , 5777
Expiration Date	12 months from date of receipt.