

Product Datasheet

PKC theta (phospho-S695) peptide (orb339224)

Description	The peptide is used to PKC theta (pS695) antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	PRKCQ
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338830 (PKC theta (phospho-S695) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q04759
Entrez	5588
Expiration Date	12 months from date of receipt.