

Product Datasheet

hNRNP C (phospho-S260) peptide (orb339225)

Description	The peptide is used to hNRNP C (pS260) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	HNRNPC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338831 (hNRNP C (phospho-S260) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P07910 , Q9Z204 , P17132
Entrez	15381 , 3183
Expiration Date	12 months from date of receipt.