

Product Datasheet

CDC37 (phospho-S13) peptide (orb339226)

Description	The peptide is used to CDC37 (pS13) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CDC37
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338832 (CDC37 (phospho-S13) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16543 , Q61081 , Q63692
Entrez	114562 , 12539 , 11140
Expiration Date	12 months from date of receipt.