

Product Datasheet

Lupus La (phospho-S366) peptide (orb339227)

Description	The peptide is used to Lupus La (pS366) antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	SSB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338833 (Lupus La (phospho-S366) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05455
Entrez	6741
Expiration Date	12 months from date of receipt.