

Product Datasheet

BRCA1 (phospho-S1423) peptide (orb339241)

Description	The peptide is used to BRCA1 (pS1423) antibody reactivity.
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	BRCA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338847 (BRCA1 (phospho-S1423) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P38398 , O54952
Entrez	497672 , 672
Expiration Date	12 months from date of receipt.