

Product Datasheet

JUND (phospho-S255) peptide (orb339242)

| | |
|----------------------------|---|
| Description | The peptide is used to JUND (pS255) antibody reactivity. |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | JUND |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb338848 (JUND (phospho-S255) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P15066 , P17535 , P52909 |
| Entrez | 16478 , 24518 , 3727 |
| Expiration Date | 12 months from date of receipt. |