

Product Datasheet

p38 (phospho-Y182) peptide (orb339257)

Description	The peptide is used to p38 (pY182) antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MAPK14
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338863 (p38 (phospho-Y182) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47811 , P70618 , Q16539
Entrez	1432 , 26416
Expiration Date	12 months from date of receipt.