

Product Datasheet

RB1 (phospho-S807) peptide (orb339263)

Description	The peptide is used to RB1 (pS807) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338869 (RB1 (phospho-S807) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P13405 , P33568 , P06400
Entrez	5925 , 24708 , 19645
Expiration Date	12 months from date of receipt.