

Product Datasheet

CD44 (phospho-S706) peptide (orb339265)

Description	The peptide is used to CD44 (pS706) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	CD44
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338871 (CD44 (phospho-S706) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P16070, P26051, P15379
Entrez	960, 12505, 25406
Expiration Date	12 months from date of receipt.