

Product Datasheet

CHK1 (phospho-S301) peptide (orb339266)

Description	The peptide is used to CHK1 (pS301) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CHEK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338872 (CHK1 (phospho-S301) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O35280 , O14757 , Q91ZN7
Entrez	140583 , 12649 , 1111
Expiration Date	12 months from date of receipt.