

Product Datasheet

mTOR (phospho-T2446) peptide (orb339267)

Description	The peptide is used to mTOR (pT2446) antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	MTOR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338873 (mTOR (phospho-T2446) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JLN9 , P42345 , P42346
Entrez	56717 , 56718 , 2475
Expiration Date	12 months from date of receipt.