

## Product Datasheet

### DCP1A peptide (orb339287)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | The peptide is used to DCP1A antibody reactivity.   |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Tested Applications</b> | Blocking  |
| <b>Target</b>              | DCP1A   |
| <b>Form/Appearance</b>     | Lyophilized powder  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb338893 (DCP1A antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. |
| <b>Purity</b>              | > 85%   |
| <b>Source</b>              | Synthetic   |
| <b>Uniprot ID</b>          | <a href="#">Q9NPI6</a> , <a href="#">Q91YD3</a>   |
| <b>Entrez</b>              | <a href="#">75901</a> , <a href="#">55802</a>   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |