

Product Datasheet

ZAP70 (phospho-Y319) peptide (orb339296)

Description	The peptide is used to ZAP70 (pY319) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ZAP70
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338902 (ZAP70 (phospho-Y319) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P43403 , P43404
Entrez	7535 , 22637
Expiration Date	12 months from date of receipt.