

Product Datasheet

STMN1 (phospho-S38) peptide (orb339297)

| | |
|----------------------------|---|
| Description | The peptide is used to STMN1 (pS38) antibody reactivity. |
| Reactivity | Human, Porcine, Primate, Rat |
| Tested Applications | Blocking |
| Target | STMN1 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb338903 (STMN1 (phospho-S38) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P13668 , P16949 |
| Entrez | 29332 , 3925 |
| Expiration Date | 12 months from date of receipt. |