

## Product Datasheet

### CPI17 (phospho-T38) peptide (orb339300)

<b>Description</b>	The peptide is used to CPI17 (pT38) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	PPP1R14A
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338906 (CPI17 (phospho-T38) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q99MC0</a> , <a href="#">Q91VC7</a> , <a href="#">Q96A00</a>
<b>Entrez</b>	<a href="#">114004</a> , <a href="#">68458</a> , <a href="#">94274</a>
<b>Expiration Date</b>	12 months from date of receipt.