

## Product Datasheet

### ACE2 peptide (orb339324)

<b>Description</b>	The peptide is used to ACE2 antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	ACE2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338930 (ACE2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q9BYF1</a> , <a href="#">Q8R0I0</a> , <a href="#">Q5EGZ1</a>
<b>Entrez</b>	<a href="#">70008</a> , <a href="#">59272</a> , <a href="#">302668</a>
<b>Expiration Date</b>	12 months from date of receipt.