

Product Datasheet

p53 (phospho-S46) peptide (orb339325)

Description	The peptide is used to p53 (pS46) antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	TP53
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338931 (p53 (phospho-S46) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04637
Entrez	7157
Expiration Date	12 months from date of receipt.