

## Product Datasheet

# Focal Adhesion Kinase (phospho-Y407) peptide (orb339327)

<b>Description</b>	The peptide is used to Focal Adhesion Kinase (pY407) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Primate, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	PTK2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338933 (Focal Adhesion Kinase (phospho-Y407) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P34152</a> , <a href="#">O35346</a> , <a href="#">Q05397</a>
<b>Entrez</b>	<a href="#">5747</a> , <a href="#">14083</a> , <a href="#">25614</a>
<b>Expiration Date</b>	12 months from date of receipt.