

Product Datasheet

AKT2 (phospho-S474) peptide (orb339339)

Description	The peptide is used to AKT2 (pS474) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	AKT2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338945 (AKT2 (phospho-S474) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47197 , P31751 , Q60823
Entrez	11652 , 208 , 25233
Expiration Date	12 months from date of receipt.