

Product Datasheet

IKB alpha (phospho-Y42) peptide (orb339340)

Description	The peptide is used to IKB alpha (pY42) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Tested Applications	Blocking
Target	NFKBIA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338946 (IKB alpha (phospho-Y42) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z1E3 , Q63746 , P25963
Entrez	25493 , 4792 , 18035
Expiration Date	12 months from date of receipt.