

## Product Datasheet

### SHCA (phospho-Y427) peptide (orb339349)

<b>Description</b>	The peptide is used to SHCA (pY427) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	SHC1
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338955 (SHCA (phospho-Y427) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P98083</a> , <a href="#">Q5M824</a> , <a href="#">P29353</a>
<b>Entrez</b>	<a href="#">6464</a> , <a href="#">20416</a> , <a href="#">85385</a>
<b>Expiration Date</b>	12 months from date of receipt.