

Product Datasheet

IRF3 (phospho-S385) peptide (orb339353)

Description	The peptide is used to IRF3 (pS385) antibody reactivity.
Reactivity	Human, Mouse, Porcine
Tested Applications	Blocking
Target	IRF3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338959 (IRF3 (phospho-S385) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q14653 , P70671
Entrez	3661 , 54131
Expiration Date	12 months from date of receipt.