

Product Datasheet

XIAP (phospho-S87) peptide (orb339355)

Description	The peptide is used to XIAP (pS87) antibody reactivity.
Reactivity	Human, Rabbit
Tested Applications	Blocking
Target	XIAP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338961 (XIAP (phospho-S87) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P98170
Entrez	331
Expiration Date	12 months from date of receipt.