

Product Datasheet

BLNK (phospho-Y96) peptide (orb339356)

Description	The peptide is used to BLNK (pY96) antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	BLNK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338962 (BLNK (phospho-Y96) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8WV28
Entrez	29760
Expiration Date	12 months from date of receipt.