

## Product Datasheet

### AKT (phospho-Y326) peptide (orb339357)

<b>Description</b>	The peptide is used to AKT (pY326) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Rat, Sheep
<b>Tested Applications</b>	Blocking
<b>Target</b>	AKT1
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb338963 (AKT (phospho-Y326) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P47196</a> , <a href="#">P31750</a> , <a href="#">P31749</a>
<b>Entrez</b>	<a href="#">207</a> , <a href="#">24185</a> , <a href="#">11651</a>
<b>Expiration Date</b>	12 months from date of receipt.