

Product Datasheet

PER2 (phospho-S662) peptide (orb339358)

Description	The peptide is used to PER2 (pS662) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	PER2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338964 (PER2 (phospho-S662) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O54943 , O15055
Entrez	18627 , 8864
Expiration Date	12 months from date of receipt.