

Product Datasheet

MKK1 (phospho-T386) peptide (orb339360)

Description	The peptide is used to MKK1 (pT386) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	MAP2K1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338966 (MKK1 (phospho-T386) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02750 , P31938 , Q01986
Entrez	26395 , 170851 , 5604
Expiration Date	12 months from date of receipt.