

Product Datasheet

MSK1 (phospho-S212) peptide (orb339361)

Description	The peptide is used to MSK1 (pS212) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	RPS6KA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338967 (MSK1 (phospho-S212) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8C050 , O75582
Entrez	9252 , 73086
Expiration Date	12 months from date of receipt.