

Product Datasheet

Inhibin beta E peptide (orb339540)

Description	The peptide is used to Inhibin beta E antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	INHBE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339146 (Inhibin beta E antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08717 , P58166 , O88959
Entrez	16326 , 83711 , 83729
Expiration Date	12 months from date of receipt.