

Product Datasheet

LHFPL1 antibody (FITC) (orb33955)

| | |
|------------------------|--|
| Description | LHFPL1 antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | FITC |
| Immunogen | Recombinant Human LHFPL tetraspan subfamily member 1 protein (21-85AA) |
| Target | LHFPL1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q86WI0 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |