

Product Datasheet

NPY5R peptide (orb339570)

Description	The peptide is used to NPY5R antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	NPY5R
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339176 (NPY5R antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q15761
Entrez	4889
Expiration Date	12 months from date of receipt.