

Product Datasheet

EGFR (phospho-Y1197) peptide (orb339585)

Description	The peptide is used to EGFR (pY1197) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	EGFR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339191 (EGFR (phospho-Y1197) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q01279 , P00533
Entrez	1956 , 13649
Expiration Date	12 months from date of receipt.