

## Product Datasheet

### Human GPNMB ELISA Kit (orb339662)

<b>Description</b>	Human Transmembrane glycoprotein NMB (GPNMB) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.78 ng/mL-50 ng/mL
<b>Target</b>	GPNMB
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q14956</a>
<b>Sensitivity</b>	0.195 ng/mL
<b>Expiration Date</b>	Please enquire.