

Product Datasheet

Human ANGPT1 ELISA Kit (orb339703)

Description	Human Angiopoietin-1 (ANGPT1) ELISA kit
Reactivity	Human
Range	1.25 ng/mL-80 ng/mL
Target	ANGPT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	Q15389
Sensitivity	0.091 ng/mL
Expiration Date	Please enquire.