

## Product Datasheet

### Human PDGF-BB ELISA Kit (orb339719)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human Platelet-Derived Growth Factor-BB, PDGF-BB ELISA kit |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                                     |
| <b>Target</b>            | PDGF-BB  |
| <b>Note</b>              | For research use only                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                               |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                          |
| <b>Sensitivity</b>       | 7.81 pg/mL   |
| <b>Expiration Date</b>   | Please enquire.  |