

## Product Datasheet

### Human VEGFR-2/Flk-1 ELISA Kit (orb339730)

<b>Description</b>	Human Vascular endothelial cell growth factor receptor 2, VEGFR-2/Flk-1 ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.39 ng/mL-25 ng/mL
<b>Target</b>	VEGFR-2/Flk-1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P35968</b>
<b>Sensitivity</b>	0.098 ng/mL
<b>Expiration Date</b>	Please enquire.