

## Product Datasheet

### Mouse ADP ELISA Kit (orb339732)

<b>Description</b>	Mouse adiponectin, ADP ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156 ng/ml-10 ng/ml.
<b>Target</b>	ADP
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, cell lysates.
<b>Uniprot ID</b>	<a href="#">Q60994</a>
<b>Sensitivity</b>	0.039 ng/ml
<b>Expiration Date</b>	Please enquire.