

Product Datasheet

Mouse ANGPTL3 ELISA Kit (orb339734)

Description Mouse Angiopoietin-related protein 3 (ANGPTL3) ELISA kit

Reactivity Mouse

Range 31.25 pg/mL-2000 pg/mL

Target ANGPTL3

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID [Q9R182](#)

Sensitivity 10.24 pg/mL

Expiration Date Please enquire.