

## Product Datasheet

### Mouse DKK1 ELISA Kit (orb339739)

<b>Description</b>	Mouse Dickkopf-related protein 1 (DKK1) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	78.125 pg/mL-5000 pg/mL
<b>Target</b>	DKK1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>O54908</b>
<b>Sensitivity</b>	13.127 pg/mL
<b>Expiration Date</b>	Please enquire.