

## Product Datasheet

### Human ALCAM ELISA Kit (orb339780)

<b>Description</b>	Human Activated Leukocyte Cell Adhesion Molecule, ALCAM ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.039 ng/mL-2.5 ng/mL
<b>Target</b>	ALCAM
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q13740</a>
<b>Sensitivity</b>	0.0097 ng/mL
<b>Expiration Date</b>	Please enquire.