

Product Datasheet

Human CCL5 ELISA Kit (orb339788)

Description	Human C-C motif chemokine 5 (CCL5/D17S136E/SCYA5) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CCL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P13501
Sensitivity	2.307 pg/mL
Expiration Date	Please enquire.