

Product Datasheet

Human CXCL16 ELISA Kit (orb339790)

| | |
|--------------------------|-------------------------------------------------|
| Description | Human CXC-chemokine ligand 16, CXCL16 ELISA Kit |
| Reactivity | Human |
| Range | 46.875 pg/mL-3000 pg/mL |
| Target | CXCL16 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q9H2A7 |
| Sensitivity | 3.603 pg/mL |
| Expiration Date | Please enquire. |